

Specifications	Mi151
Range pH	-2.00 to 16.00 pH
mV	±699.9 mV / ±1999 mV
Temp	-20.0 to 120.0°C / -4.0 to 248.0°F
Resolution pH	0.01 pH
mV	0.1 mV / 1 mV
Temp	0.1°C / 0.1°F
Accuracy pH	±0.01 pH
(@20°C) mV	±0.2 mV / 1 mV
Temp	±0.4°C / ±0.8°F
Typical EMC pH	±0.02 pH
Deviation mV	±0.2 mV / 1 mV
Temp	±0.4°C / ±0.8°F
pH Automatic Calibration	1 or 2 point calibration with 7 memorized buffers
Offset Calibration	±1 pH
Slope Calibration	from 80 to 108%
Temperature Compensation	automatic from -20.0 to 120.0°C / -4.0 to 248.0°F or manual, without temperature probe
pH Electrode	MA917B/1 (included)
Temperature Probe	MA831R (included)
Environment	0 to 50°C / 32 to 122°F; max RH 95%
Input Impedance	10 ¹² Ohm
Power Supply	12 VDC power adapter (included)
Packaging dimensions	335 x 120 x 255 mm
Packaging weight	2 kg



Glass Electrode & Temperature **Probe**

Choose from our wide selection of pH and ORP electrodes at pages 6 and 50.

Custom dual level LCD

Large and easyto-read Custom dual level LCD Display with simultaneous readings and with user-friendly icons.



Accessories

MA9001 pH 1.68 buffer solution, 230 mL bottle pH 4.01 buffer solution, 230 mL bottle MA9004 MA9006 pH 6.86 buffer solution, 230 mL bottle MA9007 pH 7.01 buffer solution, 230 mL bottle MA9009 pH 9.18 buffer solution, 230 mL bottle MA9010 pH 10.01 buffer solution, 230 mL bottle MA9012 Refilling solution for double junction electrode, 230 mL bottle MA9015 Electrode storage solution, 230 mL Electrode cleaning solution, 230 mL MA9016

Temperature probe

MA9112 pH 12.45 buffer solution, 230 mL bottle MA9310 12 VDC Adapter, 220 V

MA9311 12 VDC Adapter, 110 V MA9315 Electrode Holder

MA917B/1 Glass body, double junction refillable pH electrode

MA924B/1 ±2000 mV Glass ORP electrode, refillable with BNC connector and 1

meter cable Platinum ORP electrode with 1 m SE300 cable

Ordering Information

Mi151 is supplied complete with:

- MA917B/1 Double junction refillable pH electrode
- MA831R Temperature Probe
- MA9315 Electrode Holder
- M10004 pH 4.01 Sachet Buffer Solution
- M10007 pH 7.01 Sachet Buffer Solution
- M10010 pH 10.01 Sachet Buffer Solution • M10016 Sachet Electrode Cleaning Solution
- MA9310 12 VDC Adapter
- Instruction manual

Contact:

MA831R

EASTERN ENERGY CO.,LTD. (HEAD OFFICE)

40/4 Vitoondumri Rd., Banbueng, Banbueng, Chonburi 20170 Tel: 0-3844-6117 Fax: 0-3844-6200 Email: info@ete.co.th

www.ete.co.th