pH/ORP

3

MEM

2

CE

Mi150 pH/Temperature Laboratory Bench Meter

Mi150 is an advanced pH/Temp microprocessor-based bench meter. It is ideal for students and technicians who need fast and reliable measurements.

This meter is provided with a series of new diagnostic features which add an entirely new dimension to the measurement of pH, by allowing the user to dramatically improve the reliability of the measurement:

- Automatic Temperature Compensation (ATC) for good accuracy under fluctuating temperatures;
- Easy to read large custom LCD;
- Easy and Quick Push-button Calibration
- 7 memorized buffers (pH 1.68, 4.01, 6.86, 7.01, 9.18, 10.01 and 12.45) for calibration;
- Messages on the LCD to make the calibration easy and accurate;
- User-selectable "calibration time out" to remind when a new calibration is necessary;
- Stability Indicator prompts whenever reading stabilizes.

Moreover, it offers an extended temperature range from -20°C (-4°F) to 120°C (248°F), using the MA831R interchangeable temperature probe.

Specifications	Mi150
· · ·	
Range pH	-2.00 to 16.00 pH
Temp	-20.0 to 120.0°C / -4.0 to 248.0°F
Resolution pH	0.01 pH
Temp	0.1°C / 0.1°F
Accuracy pH	±0.01 pH
(@20°C / 68°F) Temp	±0.4°C / ±0.8°F
Typical EMC pH	±0.02 pH
Deviation Temp	±0.4°C / ±0.8°F
pH Automatic Calibration	1 or 2 point calibration with 7 memorized buffers
Offset Calibration	±1 pH
Slope Calibration	from 80 to 108%
Temperature	automatic from -20.0 to 120.0°C / -4.0 to 248.0°F
Compensation	or manual, without temperature probe
pH Electrode	MA917B/1 (included)
Temperature Probe	MA831R (included)
Environment	0 to 50°C / 32 to 122°F; max RH 95%
Input Impedance	10 ¹² Ohm
Power Supply	12 VDC power adapter (included)
Packaging dimensions	335 x 120 x 255 mm
Packaging weight	2 kg

Electrode cleaning solution,

12 VDC Adapter, 220 V

12 VDC Adapter, 110 V

Electrode Holder

Temperature probe

pH electrode

pH 12.45 buffer solution, 230 mL



230 mL bottle

. bottle

pH 1.68 buffer solution, 230 mL bottle MA9016 MA9004 pH 4.01 buffer solution, 230 mL bottle MA9006 pH 6.86 buffer solution, 230 mL bottle MA9112 MA9007 pH 7.01 buffer solution, 230 mL bottle MA9009 pH 9.18 buffer solution, 230 mL bottle MA9310 MA9010 pH 10.01 buffer solution, 230 mL bottle MA9311 MA9012 Refilling solution for double junction MA9315 electrode, 230 mL bottle MA917B/1 Glass body, double junction refillable MA9015 Electrode storage solution, 230 mL bottle MA831R

Contact :

Accessories MA9001

EASTERN ENERGY CO., LTD. (HEAD OFFICE)

40/4 Vitoondumri Rd., Banbueng, Banbueng, Chonburi 20170 Tel: 0-3844-6117 Fax: 0-3844-6200 Email: info@ete.co.th

Glass Electrode & Temperature Probe

Choose from our wide selection of pH and ORP electrodes at pages 6 and 50.

Innovative Design

Compact-size ergonomic design with electrode holder that can hold multiple electrodes & probes.



Ordering Information

Mi150 is supplied complete with:

- MA917B/1 Double junction refillable pH electrode
- MA831R Temperature Probe
- MA9315 Electrode Holder
- M10004 pH 4.01 Sachet Buffer Solution
- M10007 pH 7.01 Sachet Buffer Solution
- M10010 pH 10.01 Sachet Buffer Solution
- M10016 Sachet Electrode Cleaning Solution MA9310 12 VDC Adapter
- Instruction manual

www.ete.co.th