## **USB LOGGER w/PDF REPORT**







Ordering Code VZ88181AZ, 88181 logger (\*\* Note, 10pcs loggers included)

## Single Use Cold Chain Temperature Logger w/PDF report 88181

- Plug & Play. Software and USB driver are not required
- Automatically generate PDF report. Report contains summary, graph and detailed data
- Completely user re-programmable: 5 languages, sampling interval, start delay, alarm delay, unit, alarm range, alarm type, password and user name
- 4 LEDs provide instant indicator : REC, Low Bat., High /Low temperature alarm
- File encrypted technology to secure data integrity



















PC Link	USB 2.0	IP65	LED Alarm	
	Model			88181
Temp. range		N <mark>TC</mark> therm <mark>ist</mark> or, -30~70°C , -22~158°F		
Temp. resolution		0.1°C , 0.1°F		
Temp. Accuracy				+/-0.5°C
Logging Type				Single Use
Pouched size				105(L)*55(W)*2.5(T) mm
Unpouched size				80(L)*43(W)*2.5(T) mm
Sampl	ing points			8192 readings
Shelf L	.ife			12 months shelf time
Operat	ting temp.	-30~	70°C (Log	ging status); room temp.(PC status)
Operating RH%		Humidity < 80%		
Storage temp.		-40~85°C		
Storage RH%		Humidity < 90%		
Weight		~10g		
Battery		1PC 3.0V CR2032 (installed before shipment)		
Sampling interval		30 seconds, 5,10,30,60,90,120 minutes		
Start d	elay		C	), 5,30,45,60, 90,120 minutes
Alarm	range		-20~ -10°C, -	.10~0°C, 2~8°C, 0~15°C, 0~25°C, 15~25°C, others
Alarm delay		0, 5,30,45,60, 90,120 minutes		
Alarm type		Single Event, Cumulative, Disable		
Operation keys			2 Key	rs (Start/Stop key & Mark key)
LED indicator			REC	, Low Bat., High alarm, Low alarm

## Application:

Standard package

During transportation, record time and temperature for vaccines, biologics, pharmaceuticals and other perishable & temperature sensitive products

10pcs logger in one box, manual, paper box



